

ERGOSTAT®
PREP & PACK TABLES



 **SKYTRON**®
TRUSTED PARTNER | PROVEN SOLUTIONS



Ergonomic stand or sit flexibility for your sterile processing needs.

ErgoStat® Prep & Pack tables are designed to alleviate the physical stress placed on technicians performing the inspection, assembly and packaging in your Central Sterile Supply Department. Standing in one position for hours on end, or using a workstation designed for someone else, can lead to uncomfortable neck or back pain and unsatisfied employees. Put your facility's well-being in mind, and choose a Prep & Pack table that easily adjusts to the various physical attributes and requirements of the human body – ultimately achieving better health and employee longevity.



STAND
OR
SIT

NON-POROUS
WORK
SURFACE

EASY
ACCESSORY
ATTACHMENT

DIFFERENTIATING FEATURES

- Non-porous solid work surface
- 26 ½" to 41 ½" height adjustable work surface
- 43" vertical riser frame raises and lowers with the base
- Powered and manual (pin-adjustable) models available
- Ergonomically correct alignment of computer monitors, lighting, bins and pegboards
- Heavy duty inner/outer leg extrusions for increased stability
- 750 lb dynamic and 1,000 lb static weight capacity



OPTIONAL ACCESSORIES



Akro Bin



Bin Rail,
30" or 36"



Cantilever
Light Bracket



Vertical CPU
Holder



Electrical Strip,
30" or 36"



Waste Basket
and Holder



Footrest, 60" or 72"



Flat Screen
Keyboard Mount



Flat Screen Arm



Keyboard and
Mouse Tray



High Performance
Task Light,
60" or 72"



Lower Shelf,
18" X 60"
or 18" X 72"



Magnetic Marker
Board, 18" X 30"
or 18" X 36"



Multifunction
Shelf, 15" X 30"
or 15" X 36"




Magnification Light



Hanging Pedestal,
6" Box or 6" Box
and 12" File



Riser Frame



Spool Holder,
30" or 36"

ErgoStat® Prep & Pack Tables



SPECIFICATIONS	ERGOSTAT 36 X 72 E	ERGOSTAT 36 X 72 M	ERGOSTAT 36 X 60 E	ERGOSTAT 36 X 60 M	ERGOSTAT 30 X 72 E	ERGOSTAT 30 X 72 M	ERGOSTAT 30 X 60 E	ERGOSTAT 30 X 60 M
Frame Height	81"	81"	81"	81"	81"	81"	81"	81"
Width	73½"	73½"	61½"	61½"	73½"	73½"	61½"	61½"
Depth	39"	39"	39"	39"	33"	33"	33"	33"
Desk Height	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"
Height Adjustment	Power	Manual	Power	Manual	Power	Manual	Power	Manual

Printed with agra-based inks, water based coatings on FSC®-certified paper.



TRUSTED PARTNER | PROVEN SOLUTIONS



5085 Corporate Exchange Blvd. SE
 Grand Rapids, MI 49512
 1.800.SKYTRON (759-8766)
 email: info@skytron.us
www.skytron.us

About Skytron Skytron is a privately held company specializing in capital equipment that drives efficiency in today's healthcare facilities. Our surgical tables, surgical & exam lighting, booms, video integration systems and other products help hospitals around the world work better. To learn more, visit our website at www.skytron.us.