



1602 ESSENTIA

SPECIFICATIONS

Length w/o head rest	68 1/2"	1,740mm
Length with head rest	78 3/4"	2,000mm
Width (back / seat plate)	19 5/8"	500mm
Width with side rails	21 7/8"	555mm

Adjustments via pendant control:

Height	27 3/16" to 40 15/16"	690-1,040mm
Lateral tilt (left/right)	23°	
Trendelenburg	28°	
Reverse Trendelenburg	28°	
Back up/down	90°/30°	
Leg down	100°	
Return to level (anatomically correct return to horizontal position from articulated positions)		
Flex/Reflex	30°/90°	
Base braking (lock/unlock)		

Adjustments via manual control:

Top Rotation	180°
Head Section	60° up/90° down

Max. permissible load:	500 lb (227kg)	Lift
	500 lb (227kg)	Articulation

Battery powered surgical table features 180° top rotation and 500 lb capacity. The table is equipped with a removable 15' power cord and four conductive casters. Safety features include a four-point braking system with push button operation and manual override. Battery power ensures table operation in case of power loss. Handheld pendant permits control at any point around the table's perimeter.



5085 Corporate Exchange Blvd. SE
 Grand Rapids, MI 49512
 1.800.SKYTRON (759-8766)
 email: info@skytron.us
www.skytron.us

About Skytron Skytron is a privately held company specializing in capital equipment that drives efficiency in today's healthcare facilities. Our surgical tables, surgical & exam lighting, booms, video integration systems and other products help hospitals around the world work better. To learn more, visit our website at www.skytron.us.