






Refer to SKYTRON Mounting Structure Pre-Installation Guide for appropriate mounting structure construction and a definition of SKYTRON's and Customer's pre-installation responsibilities.  
 Final model configuration requires SKYTRON technical review and approval to determine carrier configuration and equipment load capacities.

FIXTURE WEIGHT: 240 LBS. MOMENT LOAD: 1100 FT LBS.

| Ceiling height requirements |   |   |
|-----------------------------|---|---|
|                             | Maximum   | Minimum   |
| Skylink Fixed Arm           | 120" [3048mm] =<br>77" [1956mm] bottom of SkyLink | 106" [2692mm] =<br>77" [1956mm] bottom of SkyLink |
| Monitor Spring Arm          | 120" [3048mm] =<br>80" [2032mm] center of monitor | 115" [2921mm] =<br>78" [1981mm] center of monitor |
| Skylink Spring Arm          | 120" [3048mm] =<br>79" [2007mm] bottom of SkyLink | 111" [2819mm] =<br>79" [2007mm] bottom of SkyLink |



| DATE      | REV | REVISION HISTORY   |
|-----------|-----|--|
| 6/20/2016 | 2   | Fixed typo to 9-1/2" measurement, fixed wrong number from 1866mm to 1888mm |