


FIXTURE WEIGHT: 813 LBS.    MOMENT LOAD: 4647 FT. LBS.


Ceiling height requirements		
	Maximum	Minimum
Monitor	120" [3048mm] = 77" [1956mm] center of monitor	113" [2870mm] = 71" [1803mm] center of monitor
Skylink	120" [3048mm] = 79" [2007mm] bottom of SkyLink	113" [2870mm] = 79" [2007mm] bottom of SkyLink

TEC-H-0066 REV 1  
12/15/2015



FIXTURE WEIGHT: 933 LBS.    MOMENT LOAD: 5370 FT. LBS.


Ceiling height requirements		
	Maximum	Minimum
Monitor	120" [3048mm] = 77" [1956mm] center of monitor 80" [2032mm] BOM	113" [2870mm] = 70" [1778mm] center of monitor 80" [2032mm] BOM
Skylink	120" [3048mm] = 79" [2007mm] bottom of SkyLink 80" [2032mm] BOM	113" [2870mm] = 75" [1905mm] bottom of SkyLink 80" [2032mm] BOM





Refer to SKYTRON Mounting Structure Pre-Installation Guide for appropriate mounting structure construction and a definition of SKYTRON's and Customer's pre-installation responsibilities. Final model configuration requires SKYTRON technical review and approval to determine carrier configurations and equipment load capacities.


FIXTURE WEIGHT: 1031 LBS. MOMENT LOAD: 5754 FT. LBS.


Ceiling height requirements		
	Maximum	Minimum
Monitor	120" [3048mm] = 77" [1956mm] center of monitor 80" [2032mm] BOM	118" [2997mm] = 75" [1905mm] center of monitor 80" [2032mm] BOM
Skylink	120" [3048mm] = 79" [2007mm] bottom of SkyLink 80" [2032mm] BOM	118" [2997mm] = 80" [2032mm] bottom of SkyLink 80" [2032mm] BOM





**3EC-L**  61" [1550mm]


**3EC**  51" [1300mm]


**32F-L**  FIXED HEIGHT 39 1/2" [1000mm]


**32F**  FIXED HEIGHT 31 1/2" [800mm]


**32F-C**  FIXED HEIGHT 23 1/2" [600mm]


**32S-H**  MANUAL HEIGHT ADJUSTABLE 31 1/2" [800mm]


**32S**  MANUAL HEIGHT ADJUSTABLE 31 1/2" [800mm]


**32E-XL**  POWER HEIGHT ADJUSTABLE 39 1/2" [1000mm]

**32E**  POWER HEIGHT ADJUSTABLE 31 1/2" [800mm]


**H2-300-03-A ECT Adapter**  51" [1300mm] **3T2**

**3T2-C**  41 1/2" [1050mm]

**3LF**  FIXED HEIGHT 39 1/2" [1000mm]

**3AFC-L**  
**3AFC-M**  
**3AFC-S**  
**3AFC-H**  MANUAL HEIGHT ADJUSTABLE 36" [910mm]

NO LIGHT  
REQUIRES CAP  
H2-300-04-A



**MUST ADD SPINDLE TO SUPPORT LIGHTS**

PART NO.	DESCRIPTION
H2-201-44	Light Spindle, 3ECT, NL25
H2-201-45	Light Spindle, 3ECT, NL2525
H2-201-42	Light Spindle, 3ECT, AUA Single [TV or NON-TV]
H2-201-43	Light Spindle, 3ECT, AUA Dual [TV or NON-TV]