



FIXTURE WEIGHT: 1355 LBS. MOMENT LOAD: 6846 FT. LBS.

Dimensions are based on lower arm carrier/UB

Ceiling height requirements		
	Maximum	Minimum
Carrier Fixed arm	120" [3048mm] = 71" [1804mm] top of carrier	100" [2540mm] = 70" [1778mm] top of carrier
UB Spring and elevato arm	120" [3048mm] = 78" [2032mm] bottom of UB	120" [3048mm] = 78" [1981mm] bottom of UB
Carrier Spring and elevato arm	120" [3048mm] = 71" [1804mm] top of carrier	102" [2591mm] = 70" [1778mm] top of carrier



